



Product Datasheet

CDCA3 antibody (orb223253)

Description	Rabbit polyclonal antibody to CDCA3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q99618
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using CDCA3 antibody.



Immunohistochemistry analysis of human tissue using CDCA3 antibody.