



Product Datasheet

COX17 antibody (orb223261)


Description	Rabbit polyclonal antibody to COX17
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q14061
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of 3T3 cells (Lane...



Western blot analysis of HeLa cells HepG...



Immunohistochemistry analysis of human l...