



## Product Datasheet

### ELOA1 antibody (orb223281)

|                            |                                                    |
|----------------------------|----------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ELOA1                |
| <b>Species/Host</b>        | Rabbit                                             |
| <b>Reactivity</b>          | Human, Mouse, Rat                                  |
| <b>Conjugation</b>         | Unconjugated                                       |
| <b>Tested Applications</b> | IHC, WB                                            |
| <b>Immunogen</b>           | Internal                                           |
| <b>Concentration</b>       | 1mg/ml                                             |
| <b>Note</b>                | For research use only                              |
| <b>Clonality</b>           | Polyclonal                                         |
| <b>Uniprot ID</b>          | <a href="#">Q14241</a>                             |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:1000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                    |



Immunohistochemistry analysis of human b...



Western blot analysis of Jurkat cells an...