



Product Datasheet

HOXA11 antibody (orb223301)

Description	Rabbit polyclonal antibody to HOXA11
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P31270, P31277
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HUVEC cel...



Western blot analysis of RAW264.7 cells ...