



Product Datasheet

STK32C antibody (orb223312)

Description	Rabbit polyclonal antibody to STK32C also known as PKE and Yet another novel kinase 3 whose expressi...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q86UX6
Dilution Range	WB: 1:500-3000, IHC-P: 1:50-100, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells us...



Immunohistochemistry analysis of human b...