


Product Datasheet

RAD54B antibody (orb223321)

Description	Rabbit polyclonal antibody to RAD54B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9Y620
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Immunohistochemistry analysis of human b...