



Product Datasheet

TACC3 antibody (orb223324)


Description	Rabbit polyclonal antibody to TACC3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9Y6A5
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cells...



Western blot analysis of 293 cells and H...



Immunohistochemistry analysis of human I...