



Product Datasheet

HTR5B antibody (orb223331)

Description	Rabbit polyclonal antibody to HTR5B
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P35365
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HUVEC cells...



Western blot analysis of COLO cells using...