


Product Datasheet

HTR7 antibody (orb223332)

Description	Rabbit polyclonal antibody to HTR7
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P34969
Dilution Range	IF/ICC: 1:100-1:500, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of COS-7
cel...