


Product Datasheet

CELSR3 antibody (orb223339)

Description	Rabbit polyclonal antibody to CELSR3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9NYQ7
Dilution Range	IF/ICC: 1:100-1:500, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of A549
cell...