



Product Datasheet

CXCR7 antibody (orb223343)

Description	Rabbit polyclonal antibody to CXCR7
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P25106
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of COS-7 cells...



Western blot analysis of COS-7 cells using CXCR7 antibody.