



Product Datasheet

CHRM2 antibody (orb223383)

Description	Rabbit polyclonal antibody to CHRM2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P08172
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of LOVO cell...



Western blot analysis of MCF-7 cells and...