



Product Datasheet

BAI1 antibody (orb223387)


Description	Rabbit polyclonal antibody to BAI1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O14514
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells us...



Immunofluorescence analysis of MCF-7 cel...



Immunohistochemistry analysis of human b...