



Product Datasheet

GPR137C antibody (orb223407)


Description	Rabbit polyclonal antibody to GPR137C
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8N3F9
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cells...



Western blot analysis of HepG2 cells using GPR137C antibody...



Immunohistochemistry analysis of human brain tissue...