



## Product Datasheet

### GHRHR antibody (orb223411)

|                            |                                                      |
|----------------------------|------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to GHRHR                  |
| <b>Species/Host</b>        | Rabbit                                               |
| <b>Reactivity</b>          | Human, Mouse, Rat                                    |
| <b>Conjugation</b>         | Unconjugated                                         |
| <b>Tested Applications</b> | IF, WB                                               |
| <b>Immunogen</b>           | C-terminal                                           |
| <b>Concentration</b>       | 1mg/ml                                               |
| <b>Note</b>                | For research use only                                |
| <b>Clonality</b>           | Polyclonal                                           |
| <b>Uniprot ID</b>          | <a href="#">Q02643</a>                               |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:1000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                      |



Western blot analysis of LOVO cells using...



Immunofluorescence analysis of HUVEC cells...