



Product Datasheet

GPR108 antibody (orb223416)

Description	Rabbit polyclonal antibody to GPR108
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9NPR9
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 cells using GPR108 antibody.



Immunofluorescence analysis of HUVEC cells using GPR108 antibody.