



Product Datasheet

GPRC5B antibody (orb223432)

Description	Rabbit polyclonal antibody to GPRC5B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9NZH0
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HUVEC cells 293...



Immunofluorescence analysis of MCF-7 cells...