


Product Datasheet

LGR4 antibody (orb223448)

Description	Rabbit polyclonal antibody to LGR4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9BXB1
Dilution Range	WB: 1:500-1:3000, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells us...