


Product Datasheet

LPHN1 antibody (orb223450)

Description	Rabbit polyclonal antibody to LPHN1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O94910
Dilution Range	IF/ICC: 1:100-1:500, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of LOVO
cell...