


Product Datasheet

MTNR1A antibody (orb223459)

Description	Rabbit polyclonal antibody to MTNR1A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P48039
Dilution Range	IF/ICC: 1:100-1:500, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of HepG2
cel...