


Product Datasheet

GPROPDR antibody (orb223462)

Description	Rabbit polyclonal antibody to GPROPDR
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q7PUA8
Dilution Range	WB: 1:500-1:3000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using GPROPDR antibody.