


Product Datasheet

OR10C1 antibody (orb223469)

Description	Rabbit polyclonal antibody to OR10C1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q5SUN7 , Q96KK4
Dilution Range	WB: 1:500-1:3000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of LOVO cells using OR10C1 antibody.