



Product Datasheet

OR2T11 antibody (orb223475)

Description	Rabbit polyclonal antibody to OR2T11
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8NH01
Dilution Range	WB: 1:500-1:3000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of LOVO cell...



Western blot analysis of MCF-7 cells and...