


Product Datasheet

OR5AP2 antibody (orb223492)

Description	Rabbit polyclonal antibody to OR5AP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8NGF4
Dilution Range	IF/ICC: 1:100-1:500, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of HeLa
cell...