


Product Datasheet

OR5AS1 antibody (orb223493)

Description	Rabbit polyclonal antibody to OR5AS1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8N127
Dilution Range	IF/ICC: 1:100-1:500, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of LOVO
cell...