


Product Datasheet

OR6C68 antibody (orb223494)

Description	Rabbit polyclonal antibody to OR6C68
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	A6NDL8
Dilution Range	WB: 1:500-1:3000, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HUVEC cells using OR6C68 antibody.