



Product Datasheet

OR10D4 antibody (orb223504)

Description	Rabbit polyclonal antibody to OR10D4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8NGN7
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO cells using OR10D4 antibody.



Immunofluorescence analysis of LOVO cells using OR10D4 antibody.