


Product Datasheet

OR1S1/1S2 antibody (orb223516)

Description	Rabbit polyclonal antibody to OR1S1/1S2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8NGQ3, Q8NH92
Dilution Range	WB: 1:500-1:3000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of HT-29 cells using OR1S1/1S2 antibody.