



Product Datasheet

OR2T1 antibody (orb223530)

Description	Rabbit polyclonal antibody to OR2T1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O43869
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HT-29 cells using OR2T1 antibody.



Immunofluorescence analysis of A549 cells using OR2T1 antibody.