



Product Datasheet

OR2T2/35 antibody (orb223531)

Description	Rabbit polyclonal antibody to OR2T2/35
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q6IF00, Q8NGX2
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using OR2T2/35 antibody.



Immunofluorescence analysis of A549 cells using OR2T2/35 antibody.