



## Product Datasheet

### OR4A16 antibody (orb223536)

|                            |                                                       |
|----------------------------|-------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to OR4A16                  |
| <b>Species/Host</b>        | Rabbit                                                |
| <b>Reactivity</b>          | Human                                                 |
| <b>Conjugation</b>         | Unconjugated                                          |
| <b>Tested Applications</b> | IF, WB                                                |
| <b>Immunogen</b>           | C-terminal                                            |
| <b>Concentration</b>       | 1mg/ml                                                |
| <b>Note</b>                | For research use only                                 |
| <b>Clonality</b>           | Polyclonal                                            |
| <b>Uniprot ID</b>          | <b>Q8NH70</b>                                         |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:20000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                       |



Western blot analysis of COLO cells using OR4A16 antibody.



Immunofluorescence analysis of A549 cells using OR4A16 antibody.