



Product Datasheet

OR4F4 antibody (orb223543)

Description	Rabbit polyclonal antibody to OR4F4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q96R69
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using OR4F4 antibody.



Immunofluorescence analysis of A549 cells using OR4F4 antibody.