


## Product Datasheet

### OR5P3 antibody (orb223557)

|                            |                                     |
|----------------------------|-------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to OR5P3 |
| <b>Species/Host</b>        | Rabbit                              |
| <b>Reactivity</b>          | Human                               |
| <b>Conjugation</b>         | Unconjugated                        |
| <b>Tested Applications</b> | WB                                  |
| <b>Immunogen</b>           | Internal                            |
| <b>Concentration</b>       | 1mg/ml                              |
| <b>Note</b>                | For research use only               |
| <b>Clonality</b>           | Polyclonal                          |
| <b>Uniprot ID</b>          | <b>Q8WZ94</b>                       |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, ELISA: 1:10000    |
| <b>Expiration Date</b>     | 12 months from date of receipt.     |



Western blot analysis of COLO cells using OR5P3 antibody.