


Product Datasheet

OR8H2 antibody (orb223569)

Description	Rabbit polyclonal antibody to OR8H2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8N162
Dilution Range	IF/ICC: 1:100-1:500, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of A549
cell...