


## Product Datasheet

### TAS2R8 antibody (orb223597)

|                            |                                      |
|----------------------------|--------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to TAS2R8 |
| <b>Species/Host</b>        | Rabbit                               |
| <b>Reactivity</b>          | Human                                |
| <b>Conjugation</b>         | Unconjugated                         |
| <b>Tested Applications</b> | WB                                   |
| <b>Immunogen</b>           | Internal                             |
| <b>Concentration</b>       | 1mg/ml                               |
| <b>Note</b>                | For research use only                |
| <b>Clonality</b>           | Polyclonal                           |
| <b>Uniprot ID</b>          | <b>Q9NYW2</b>                        |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, ELISA: 1:10000     |
| <b>Expiration Date</b>     | 12 months from date of receipt.      |



Western blot analysis of COLO cells using...