


Product Datasheet

OR10A5 antibody (orb223604)

Description	Rabbit polyclonal antibody to OR10A5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9H207
Dilution Range	IF/ICC: 1:100-1:500, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of MCF-7
cel...