


Product Datasheet

DLL4 antibody (orb223627)

Description	Mouse monoclonal antibody to DLL4
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	recombinant soluble human DLL4
Target	DLL4
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9NR61
Entrez	54567
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of
DLL4
transfected
HD4 cells...