


Product Datasheet

CD94 antibody (orb223631)

Description	Mouse monoclonal antibody to CD94
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Cultured human NK cells
Target	CD94
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q13241
Entrez	3824
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...