

Product Datasheet

Acetyl Lysine antibody (HRP) (orb223637)

Description	Mouse monoclonal antibody against Acetyl Lysine conjugated to HRP
Species/Host	Mouse
Reactivity	All
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Acetylated keyhole limpet hemocyanin (KLH) was used as immunogen. The cell line was screened against non acetylated KLH as a control in addition to BSA and histones
Target	Acetyl Lysine
Preservatives	Purified IgG conjugated to Horseradish Peroxidase, supplied in Phosphate buffered saline (PBS) containing 0.01% Thiomersal
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration.
Expiration Date	12 months from date of receipt.