

Product Datasheet

CD33 antibody (orb223670)

Description	Mouse monoclonal antibody to CD33
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-P
Immunogen	CD33-Fc recombinant protein
Target	CD33
Preservatives	Purified IgG supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	Cells from immunized mice were fused with the Sp2/0-Ag.14 myeloma cell line
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P20138
Entrez	945
Dilution Range	Optimal antibody dilution should be determined by titration
Expiration Date	12 months from date of receipt.