


Product Datasheet

CD33 antibody (orb223672)

Description	Mouse monoclonal antibody to CD33
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, WB
Immunogen	Human AML cells
Target	CD33
Preservatives	Purified antibody in PBS pH 7.2 plus 0.1% sodium azide
Concentration	1.56mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	p20138
Entrez	945
Dilution Range	Optimal antibody dilution should be determined by titration however as a guideline use at 2ug/ml for IHC-(fr)
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of using CD...