

Product Datasheet

CD41 antibody (FITC) (orb223674)

Description	Mouse monoclonal antibody against CD41 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC, ICC
Immunogen	Homogenized human medulloblastoma tissue.
Target	CD41
Preservatives	Purified IgG conjugated to Fluorescein using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	Cells from immunized Balb/c mice were fused with the P3X63Ag8.653 myeloma cell line
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P08514
Entrez	3689
Dilution Range	Optimal antibody dilution should be determined by titration
Expiration Date	12 months from date of receipt.