

Product Datasheet

CD45 antibody (orb223677)

Description	Mouse monoclonal antibody against CD45 conjugated to HRP
Species/Host	Mouse
Reactivity	Human
Conjugation	HRP
Tested Applications	IHC-Fr, IHC-P, IP, WB
Immunogen	Cultured T cells from an IL-2-dependent T-cell line (CA1) prepared from human peripheral blood activated with influenza virus.
Target	CD45RO
Preservatives	Purified IgG conjugated to Horseradish Peroxidase using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.01% Thiomersal
Concentration	1mg/ml
Note	For research use only
Application notes	Cells from immunized Balb/c mice were fused with the P3/NS1/1-Ag4.1 myeloma cell line
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P08575
Entrez	5788
Dilution Range	Optimal antibody dilution should be determined by titration
Expiration Date	12 months from date of receipt.