


Product Datasheet

GLP1R antibody (orb223703)

Description	Gallus monoclonal antibody to GLP1R
Species/Host	Gallus
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A synthetic peptide derived from the internal framework of glucagon-like 1 receptor
Target	GLP1R
Preservatives	Purified IgY supplied in Phosphate buffered saline (PBS) containing 0.01% Thiomersal and 1% trehalose
Concentration	1mg/ml
Note	For research use only
Isotype	IgY
Clonality	Polyclonal
Dilution Range	Optimal antibody dilution should be determined by titration.
Expiration Date	12 months from date of receipt.



Western blot analysis of using GLP1R ant...