

Product Datasheet

Lipoprotein Lipase antibody (HRP) (orb223715)

Description	Mouse monoclonal antibody against Lipoprotein Lipase conjugated to HRP
Species/Host	Mouse
Reactivity	Bovine, Gallus, Guinea pig, Human, Rat
Conjugation	HRP
Tested Applications	ELISA, FC, IF, WB
Immunogen	Purified bovine milk lipoprotein lipase protein
Target	Lipoprotein Lipase
Preservatives	Purified IgG conjugated to Horseradish Peroxidase using Innova Biosciences Lightning-Link, containing Phosphate buffered saline and 0.01% Sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P06858
Entrez	4023
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; IF: Use at 2 µg/ml WB: Use at an assay dependent dilution. Predicted molecular weight: 53 kDa. Flow Cyt: Use 1µg for 106 cells.
Expiration Date	12 months from date of receipt.