

Product Datasheet

Nestin antibody (FITC) (orb223719)

Description	Mouse monoclonal antibody against Nestin conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC, ICC
Immunogen	A 150-amino-acid fragment from the cloned human nestin conjugated to GST
Target	Nestin
Preservatives	Purified IgG conjugated to Fluorescein using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P48681
Entrez	10763
Dilution Range	Optimal antibody dilution should be determined by titration.
Expiration Date	12 months from date of receipt.