

Product Datasheet

Nestin antibody (HRP) (orb223720)

Description	Mouse monoclonal antibody against Nestin conjugated to HRP
Species/Host	Mouse
Reactivity	Human
Conjugation	HRP
Tested Applications	ICC, IHC-P, WB
Immunogen	A 150-amino-acid fragment from the cloned human nestin conjugated to GST
Target	Nestin
Preservatives	Purified IgG conjugated to Horseradish Peroxidase using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.01% Thiomersal
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P48681
Entrez	10763
Dilution Range	Optimal antibody dilution should be determined by titration. Predicted molecular weight: 177 kDa.
Expiration Date	12 months from date of receipt.