


Product Datasheet

C4d Complement fragment antibody (orb223757)

Description	Mouse monoclonal antibody to C4d Complement fragment
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P
Immunogen	Complement C4d from normal human plasma
Target	C4d Complement
Preservatives	Purified antibody containing PBS + 10mM sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	p0c0l4, p0c0l5
Entrez	100293534 , 721 , 720 , 100507685
Dilution Range	Optimal antibody dilution should be determined by titration, but as a guideline use 1:25-1:50 for IHC-P
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of using C4...