

Product Datasheet

Cathelicidin antibody (FITC) (orb223758)

Description	Mouse monoclonal antibody against Cathelicidin conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC, ICC, IP, WB
Immunogen	Partially purified human neutrophil membrane proteins (heparin Ultrogellectin affinity purification fractions in octyl glucoside)
Target	Cathelicidin [H7]
Preservatives	Purified IgG conjugated to Fluorescein using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.09% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	P3U1 myeloma;
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	p49913
Entrez	820
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try at WB: 1:1000. Predicted molecular weight: 19 kDa
Expiration Date	12 months from date of receipt.