

Product Datasheet

Complement C4d antibody (FITC) (orb223764)

Description	Mouse monoclonal antibody against Complement C4d conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	IHC-Fr
Immunogen	Complement C4d from normal human plasma
Target	C4d Complement
Preservatives	Purified IgG conjugated to Fluorescein using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	p0c0l4, p0c0l5
Entrez	100293534 , 721 , 720 , 100507685
Dilution Range	Optimal antibody dilution should be determined by titration.
Expiration Date	12 months from date of receipt.